



## P?r îndreptat p?rul - GIUBRA

### QHT1 NANOARGINT TITANIUM



✉ [bbazar@globelife.com](mailto:bbazar@globelife.com)  
☎ +39 0331 1706328  
🌐 [www.globelife.com/giubra/](http://www.globelife.com/giubra/)



### QHT1 NANOARGINT TITANIUM

- 230 ° C -. temperatura maxim?
- QHT -. Tehnologie de înc?lzire rapid? pentru realizarea rapid? ?i consistent? a temperaturii maxime
- Nano Titanium Silver : tehnologia aplicat la QHT1 plac? permite o c?lcare perfect?. Nanoparticule de argint, care sunt prezente in placi, asigura o sterilizare ?i tratament antibacterian.
- LCD : controlul digital al temperaturii de la 140 ° C pân? la 230 ° C.
- 10 set?ri diferite de temperatura : 140 ° C pân? la 230 ° C.
- Sistem plutitoare .. inclinare avansata a pl?cilor
- On / Off (Pornit / Oprit).
- Swivel Cord -. la 360 ° Roti?i cablul
- Cablu mai lung : 3 metri
- Tensiune dubl? - 110-240 Volt -. 50 Hz, 40 W
- L?ime 23 mm plac?.
- Cod de 6840.